

## 兔抗 CAP2 多克隆抗体

中文名称：兔抗 CAP2 多克隆抗体

英文名称： Anti-CAP2 rabbit polyclonal antibody

储 存： 冷冻（-20℃）

抗 原： CAP2

宿 主： Rabbit

反应种属： Human, Mouse, Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene was identified by its similarity to the gene for human adenylyl cyclase-associated protein. The function of the protein encoded by this gene is unknown. However, the protein appears to be able to interact with a denylyl cyclase-associated protein and actin.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	CAP2
<b>Immunogen:</b>	Fusion protein of human CAP2
<b>Full name:</b>	cyclase associated actin cytoskeleton regulatory protein 2
<b>SwissProt:</b>	P40123
<b>ELISA Recommended dilution:</b>	1000-5000

<b>IHC positive control:</b>	Human lung cancer and Human ovarian cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	53 kDa
<b>WB Positive control:</b>	Human fetal muscle tissue
<b>WB Recommended dilution:</b>	200-1000



