

WTAP 抗原（重组蛋白）

中文名称：WTAP 抗原（重组蛋白）

英文名称：WTAP Antigen (Recombinant Protein)

别名：Mum2

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	Wilms tumor 1 associated protein
Synonyms:	Mum2
Swissprot:	Q15007
Gene Accession:	BC000383
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The Wilms tumor suppressor gene WT1 appears to play a role in both transcriptional and posttranscriptional regulation of certain cellular genes. This gene encodes a WT1-associating protein, which is a ubiquitously expressed nuclear protein. Like WT1 protein, this protein is localized throughout the nucleoplasm as well as in speckles and partially colocalizes with splicing factors. Alternative splicing of this gene results in several transcript variants encoding three different isoforms.