

YPEL1 抗原（重组蛋白）

中文名称：YPEL1 抗原（重组蛋白）

英文名称： YPEL1 Antigen (Recombinant Protein)

别 名： FKSG3

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

Full name:	yippee like 1
Synonyms:	FKSG3
Swissprot:	O60688
Gene Accession:	BC074501
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene is located in the region associated with DiGeorge syndrome on chromosome 22. The encoded protein localizes to the centrosome and nucleolus and may play a role in the regulation of cell division.