

兔抗 DGKE 多克隆抗体

中文名称：兔抗 DGKE 多克隆抗体

英文名称： Anti-DGKE rabbit polyclonal antibody

别名： DGK; DAGK5; DAGK6; NPHS7

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： DGKE

宿主： Rabbit

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Diacylglycerol kinases are thought to be involved mainly in the regeneration of phosphatidylinositol (PI) from diacylglycerol in the PI-cycle during cell signal transduction. When expressed in mammalian cells, DGK-epsilon shows specificity for arachidonyl-containing diacylglycerol. DGK-epsilon is expressed predominantly in testis.
Applications:	WB
Name of antibody:	DGKE
Immunogen:	Synthesized peptide derived from internal of human DGKE.
Full name:	diacylglycerol kinase, epsilon 64kDa
Synonyms :	DGK; DAGK5; DAGK6; NPHS7

SwissProt:	P52429
WB Predicted band size:	64 kDa
WB Positive control:	K562 cells lysate
WB Recommended dilution:	500-3000

