

兔抗 DGKH 多克隆抗体

中文名称：兔抗 DGKH 多克隆抗体

英文名称： Anti-DGKH rabbit polyclonal antibody

别名： DGKeta

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： DGKH

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the diacylglycerol kinase (DGK) enzyme family. Members of this family are involved in regulating intracellular concentrations of diacylglycerol and phosphatidic acid. Variation in this gene has been associated with bipolar disorder. Alternatively spliced transcript variants have been identified.
Applications:	IHC, IF
Name of antibody:	DGKH
Immunogen:	Synthesized peptide derived from internal of human DGKH.
Full name:	diacylglycerol kinase, eta
Synonyms :	DGKeta

SwissProt:	Q86XP1
IHC positive control:	Human colon carcinoma tissue
IHC Recommend dilution:	50-100
IF positive control:	Hela cells
IF Recommend dilution:	100-500



