

兔抗 CAPN11 多克隆抗体

中文名称：兔抗 CAPN11 多克隆抗体

英文名称： Anti-CAPN11 rabbit polyclonal antibody

别名： calpain11

抗原： CAPN11

储存： 冷冻（-20℃） 避光

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Calpains constitute a family of intracellular calcium-dependent cysteine proteases. There are eight members in this superfamily. They consist of a variable 80 kDa subunit and an invariant 30 kDa subunit. This calpain protein appears to have protease activity and calcium-binding ability. A similar mouse protein may play a functional role in spermatogenesis and in the regulation of calcium-dependent signal transduction events during meiosis.
Applications:	WB
Name of antibody:	CAPN11

Immunogen:	Synthesized peptide derived from internal of human CAPN11.
Full name:	calpain 11
Synonyms :	calpain11
SwissProt:	Q9UMQ6
WB Predicted band size:	84 kDa
WB Positive control:	HepG2 cells, Jurkat cells and Huvec cells lysates
WB Recommended dilution:	500-3000

