

## ZMAT2 抗原（重组蛋白）

中文名称：ZMAT2 抗原（重组蛋白）

英文名称：ZMAT2 Antigen (Recombinant Protein)

别名：Snu23; Ptg-12; hSNU23

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	zinc finger, matrin-type 2
<b>Synonyms:</b>	Snu23; Ptg-12; hSNU23
<b>Swissprot:</b>	Q96NC0
<b>Gene Accession:</b>	BC056668
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The function of ZMAT2 remains unknown.