

兔抗 CAPN12 多克隆抗体

中文名称：兔抗 CAPN12 多克隆抗体

英文名称： Anti-CAPN12 rabbit polyclonal antibody

抗 原： CAPN12

储 存： 冷冻（-20℃） 避光

宿 主： Rabbit

反应种属： Human Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The calpains, calcium-activated neutral proteases, are nonlysosomal, intracellular cysteine proteases. The mammalian calpains include ubiquitous, stomach-specific, and muscle-specific proteins. The ubiquitous enzymes consist of heterodimers with distinct large, catalytic subunits associated with a common small, regulatory subunit. This gene encodes a member of the calpain large subunit family.
Applications:	WB
Name of antibody:	CAPN12
Immunogen:	Synthesized peptide derived from internal of human CAPN12.
Full name:	calpain 12

Synonyms :	0
SwissProt:	Q6ZSI9
WB Predicted band size:	81 kDa
WB Positive control:	NIH/3T3 cells, A549 cells and K562 cells lysates
WB Recommended dilution:	500-3000

