

兔抗 DGKK 多克隆抗体

中文名称：兔抗 DGKK 多克隆抗体

英文名称： Anti-DGKK rabbit polyclonal antibody

储 存： 冷冻（-20℃） 避光

抗 原： DGKK

宿 主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The protein encoded by this gene is an enzyme that phosphorylates diacylglycerol, converting it to phosphatidic acid. The encoded protein is a membrane protein and is inhibited by hydrogen peroxide. Variations in this gene have been associated with hypoadipias. [provided by RefSeq,
Applications:	IHC, IF
Name of antibody:	DGKK
Immunogen:	Synthesized peptide derived from C-terminal of human DGKK.
Full name:	diacylglycerol kinase, kappa
Synonyms :	0
SwissProt:	Q5KSL6

IHC positive control:	Human brain tissue
IHC Recommend dilution:	50-100
IF positive control:	COS-7 cells
IF Recommend dilution:	100-500



