

## 兔抗 CAPN15 多克隆抗体

中文名称: 兔抗 CAPN15 多克隆抗体

英文名称: Anti-CAPN15 rabbit polyclonal antibody

别名: SOLH

抗原: CAPN15

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein containing zinc-finger-like repeats and a calpain-like protease domain. The encoded protein may function as a transcription factor, RNA-binding protein, or in protein-protein interactions during visual system development.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	CAPN15
<b>Immunogen:</b>	Synthetic peptide of human CAPN15
<b>Full name:</b>	calpain 15
<b>Synonyms:</b>	SOLH
<b>SwissProt:</b>	O75808

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer and human lung cancer
<b>IHC Recommend dilution:</b>	40-200
<b>WB Predicted band size:</b>	117 kDa
<b>WB Positive control:</b>	Hela and A172 cell lysates
<b>WB Recommended dilution:</b>	500-2000



