

ZNF449 抗原（重组蛋白）

中文名称：ZNF449 抗原（重组蛋白）

英文名称：ZNF449 Antigen (Recombinant Protein)

别名：ZSCAN19

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	zinc finger protein 449
Synonyms:	ZSCAN19
Swissprot:	Q6P9G9
Gene Accession:	BC065938
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a nuclear protein that likely functions as a transcription factor. The protein includes an N-terminal SCAN domain, and several C2H2-type zinc finger motifs.