

ZNF483 抗原（重组蛋白）

中文名称：ZNF483 抗原（重组蛋白）

英文名称：ZNF483 Antigen (Recombinant Protein)

别名：ZSCAN48; ZKSCAN16

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	zinc finger protein 483
Synonyms:	ZSCAN48; ZKSCAN16
Swissprot:	Q8TF39
Gene Accession:	BC065738
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May be involved in transcriptional regulation.