

## 兔抗 DGKZ 多克隆抗体

中文名称：兔抗 DGKZ 多克隆抗体

英文名称：Anti-DGKZ rabbit polyclonal antibody

别名：DAGK5; DAGK6; DGK-ZETA; hDGKzeta

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DGKZ

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene belongs to the eukaryotic diacylglycerol kinase family. It may attenuate protein kinase C activity by regulating diacylglycerol levels in intracellular signaling cascade and signal transduction. Alternative splicing occurs at this locus and multiple transcript variants encoding distinct isoforms have been identified.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DGKZ
<b>Immunogen:</b>	Synthetic peptide of human DGKZ
<b>Full name:</b>	diacylglycerol kinase, zeta

<b>Synonyms :</b>	DAGK5; DAGK6; DGK-ZETA; hDGKzeta
<b>SwissProt:</b>	Q13574
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human cervical cancer and Human colon cancer
<b>IHC Recommend dilution:</b>	50-200

