

## 兔抗 DGLUCY 多克隆抗体

中文名称：兔抗 DGLUCY 多克隆抗体

英文名称：Anti-DGLUCY rabbit polyclonal antibody

别名：D-glutamate cyclase; C14orf159

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DGLUCY

反应种属：Human


标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	D-glutamate cyclase that converts D-glutamate to 5-oxo-D-proline.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	DGLUCY
<b>Immunogen:</b>	Fusion protein of human DGLUCY
<b>Full name:</b>	D-glutamate cyclase
<b>Synonyms:</b>	C14orf159
<b>SwissProt:</b>	Q7Z3D6
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	25-100

<b>WB Predicted band size:</b>	66 kDa
<b>WB Positive control:</b>	Human thyroid tissue lysate
<b>WB Recommended dilution:</b>	500-2000



**mlbio** 梅联生物  
Good elisakit producers

