

ZNF707 抗原（重组蛋白）

中文名称：ZNF707 抗原（重组蛋白）

英文名称： ZNF707 Antigen (Recombinant Protein)

别名： zinc finger protein 707

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 172-371 amino acids of human ZNF707

技术规格

Full name:	zinc finger protein 707
Swissprot:	Q96C28
Gene Accession:	BC014894
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May be involved in transcriptional regulation.