

兔抗 DHCR24 多克隆抗体

中文名称：兔抗 DHCR24 多克隆抗体

英文名称： Anti-DHCR24 rabbit ployclonal antibody

别名： DCE; SELADIN1; Nbla03646; seladin-1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： DHCR24

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a flavin adenine dinucleotide (FAD)-dependent oxidoreductase which catalyzes the reduction of the delta-24 double bond of sterol intermediates during cholesterol biosynthesis. The protein contains a leader sequence that directs it to the endoplasmic reticulum membrane. Missense mutations in this gene have been associated with desmosterolosis.
Applications:	ELISA, WB, IHC
Name of antibody:	DHCR24
Immunogen:	Fusion protein of human DHCR24

Full name:	24-dehydrocholesterol reductase
Synonyms :	DCE; SELADIN1; Nbla03646; seladin-1
SwissProt:	Q15392
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human breast cancer and human thyroid cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	60 kDa
WB Positive control:	Mouse brain tissue
WB Recommended dilution:	500-2000



