

兔抗 CAPS 多克隆抗体

中文名称：兔抗 CAPS 多克隆抗体

英文名称： Anti-CAPS rabbit polyclonal antibody

别名： CAPS1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CAPS

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a calcium-binding protein, which may play a role in the regulation of ion transport. A similar protein was first described as a potentially important regulatory protein in the dog thyroid and was termed as R2D5 antigen in rabbit. Alternative splicing of this gene generates two transcript variants.
Applications:	ELISA, WB, IHC
Name of antibody:	CAPS
Immunogen:	Synthetic peptide of human CAPS
Full name:	calcyphosine
Synonyms:	CAPS1

SwissProt:	Q13938
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human tonsil
IHC Recommend dilution:	20-100
WB Predicted band size:	21 kDa
WB Positive control:	Human pancreas tissue and Human kidney tissue lysates
WB Recommended dilution:	200-1000



