

兔抗 CAPZA2 多克隆抗体

中文名称：兔抗 CAPZA2 多克隆抗体

英文名称： Anti-CAPZA2 rabbit polyclonal antibody

别名： CAPZ; CAPP2

相关类别： 一抗

抗原： CAPZA2

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the F-actin capping protein alpha subunit family. It is the alpha subunit of the barbed-end actin binding protein Cap Z. By capping the barbed end of actin filaments, Cap Z regulates the growth of the actin filaments at the barbed end.
Applications:	ELISA, WB, IHC
Name of antibody:	CAPZA2
Immunogen:	Fusion protein of human CAPZA2
Full name:	capping protein (actin filament) muscle Z-line, alpha 2
Synonyms :	CAPZ; CAPP2

SwissProt:	P47755
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and human gastric cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	33 kDa
WB Positive control:	Hela cells, mouse brain and heart tissue
WB Recommended dilution:	200-1000



