

兔抗 CAPZA3 多克隆抗体

中文名称：兔抗 CAPZA3 多克隆抗体

英文名称：Anti-CAPZA3 rabbit polyclonal antibody

别名：Gsg3; CAPPA3

相关类别：一抗

抗原：CAPZA3

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an actin capping protein, one of the F-actin capping protein alpha subunit family. The encoded protein is predominantly localized to the neck region of ejaculated sperm, other immunohistochemical signals were found in the tail and postacrosomal regions. The encoded protein may also form heterodimers of alpha and beta subunits. This protein may be important in determining sperm architecture and male fertility.
Applications:	ELISA, IHC
Name of antibody:	CAPZA3
Immunogen:	Fusion protein of human CAPZA3

Full name:	capping protein (actin filament) muscle Z-line, alpha 3
Synonyms :	Gsg3; CAPP3
SwissProt:	Q96KX2
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	50-200

