

兔抗 CARD6 多克隆抗体

中文名称： 兔抗 CARD6 多克隆抗体

英文名称： Anti-CARD6 rabbit polyclonal antibody

别名： CINCIN1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CARD6

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a protein that contains a caspase recruitment domain (CARD), an antiparallel six-helical bundle that mediates homotypic protein-protein interactions. The encoded protein is a microtubule-associated protein that has been shown to interact with receptor-interacting protein kinases and positively modulate signal transduction pathways converging on activation of the inducible transcription factor NF- κ B.
Applications:	ELISA, IHC
Name of antibody:	CARD6
Immunogen:	Synthetic peptide of human CARD6

Full name:	caspase recruitment domain family, member 6
Synonyms :	CINCIN1
SwissProt:	Q9BX69
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human Lymphoma
IHC Recommend dilution:	25-100

