

## 兔抗 CARD17 多克隆抗体

中文名称：兔抗 CARD17 多克隆抗体

英文名称：Anti-CARD17 rabbit polyclonal antibody

别名：INCA

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CARD17

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Regulator of procaspase-1/CASP1 activation implicated in the regulation of the proteolytic maturation of pro-IL-1beta/IL1B and its release during inflammation. Inhibits the release of IL1B in response to LPS in monocytes. However, unlike CASP1, do not induce NF-kappa-B activation.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CARD17
<b>Immunogen:</b>	Synthetic peptide of human CARD17
<b>Full name:</b>	caspase recruitment domain family, member 17
<b>Synonyms :</b>	INCA

<b>SwissProt:</b>	Q5XLA6
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human liver cancer
<b>IHC Recommend dilution:</b>	50-200

