

兔抗 DHRS4 多克隆抗体

中文名称：兔抗 DHRS4 多克隆抗体

英文名称：Anti-DHRS4 rabbit polyclonal antibody

别名：EPD; PDE; ATQ1

抗原：DHRS4

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Reduces all-trans-retinal and 9-cis retinal. Can also catalyze the oxidation of all-trans-retinol with NAD P as co-factor, but with much lower efficiency. Reduces alkyl phenyl ketones and alpha-dicarbonyl compounds with aromatic rings, such as pyrimidine-4-aldehyde, 3-benzoylpyridine, 4-benzoylpyridine, menadione and 4-hexanoylpyridine. Has no activity towards aliphatic aldehydes and ketones (By similarity).
Applications:	ELISA, IHC
Name of antibody:	DHRS4

Immunogen:	Synthetic peptide corresponding to a region derived from 229-242 amino acids of human DHRS4
Full name:	dehydrogenase/reductase 4
Synonyms:	CR; NRDR; PHCR; PSCD; SDR-SRL; SDR25C1; SDR25C2; SCAD-SRL
SwissProt:	Q9BTZ2
ELISA Recommended dilution:	5000-10000
IHC Positive control:	Human ovarian cancer;
IHC Recommended dilution:	20-100

