

兔抗 CARTPT 多克隆抗体

中文名称：兔抗 CARTPT 多克隆抗体

英文名称： Anti-CARTPT rabbit polyclonal antibody

别名： CART

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CARTPT

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a secreted protein which is processed by prohormone/proprotein convertases to produce smaller , biologically active peptides. Expression of the transcript for this gene is regulated by certain drugs such as cocaine, and the encoded protein is thought to be involved in the regulation of appetite and stress. Mutations in this gene are associated with susceptibility to obesity.
Applications:	ELISA, WB, IHC
Name of antibody:	CARTPT
Immunogen:	Full length fusion protein
Full name:	CART prepropeptide

Synonyms :	CART
SwissProt:	Q16568
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	13 kDa
WB Positive control:	Rat heart tissue and Mouse brain tissue
WB Recommended dilution:	500-2000



