

兔抗 DHX8 多克隆抗体

中文名称：兔抗 DHX8 多克隆抗体

英文名称： Anti-DHX8 rabbit polyclonal antibody

别名： DEAH-box helicase 8; DDX8; Dhr2; HRH1; PRP22; PRPF22

抗原： DHX8

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

IHC Recommend dilution:	100-300
IHC positive control:	Human esophagus cancer and Human gastric cancer
ELISA Recommended dilution:	5000-10000
SwissProt:	Q14562
Synonyms:	DDX8; Dhr2; HRH1; PRP22; PRPF22
Full name:	DEAH-box helicase 8
Immunogen:	Fusion protein of human DHX8
Name of antibody:	DHX8
Applications:	ELISA, IHC
Background:	This gene is a member of the DEAH box polypeptide

de family. The encoded protein contains the DEAH (Asp-Glu-Ala-His) motif which is characteristic of all DEAH box proteins, and is thought to function as an ATP-dependent RNA helicase that regulates the release of spliced mRNAs from spliceosomes prior to their export from the nucleus. This protein may be required for the replication of human immunodeficiency virus type 1 (HIV-1). Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2014]



