

## 兔抗 CASP8AP2 多克隆抗体

中文名称：兔抗 CASP8AP2 多克隆抗体

英文名称：Anti-CASP8AP2 rabbit polyclonal antibody

别名：CED-4; FLASH; RIP25

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CASP8AP2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This protein is highly similar to FLASH, a mouse apoptotic protein identified by its interaction with the death-effector domain (DED) of caspase 8. Studies of FLASH protein suggested that this protein may be a component of the death-inducing signaling complex that includes Fas receptor, Fas-binding adapter FADD, and caspase 8, and plays a regulatory role in Fas-mediated apoptosis. Alternative splicing results in multiple transcript variants encoding the same protein.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CASP8AP2

<b>Immunogen:</b>	Synthetic peptide of human CASP8AP2
<b>Full name:</b>	caspase 8 associated protein 2
<b>Synonyms :</b>	CED-4; FLASH; RIP25
<b>SwissProt:</b>	Q9UKL3
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human brain
<b>IHC Recommend dilution:</b>	25-100

