

兔抗 CASP10 多克隆抗体

中文名称：兔抗 CASP10 多克隆抗体

英文名称：Anti-CASP10 rabbit polyclonal antibody

别名：caspase 10 MCH4; ALPS2; FLICE2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CASP10

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a protein which is a member of the cysteine-aspartic acid protease (caspase) family. Sequential activation of caspases plays a central role in the execution-phase of cell apoptosis. Caspases exist as inactive proenzymes which undergo proteolytic processing at conserved aspartic residues to produce two subunits, large and small, that dimerize to form the active enzyme. This protein cleaves and activates caspases 3 and 7, and the protein itself is processed by caspase 8. Mutations in this gene are associated with type IIA autoimmune lymphoproliferative syndrome, non-Hodgkin lymphoma and gastric cancer. Altern

	actively spliced transcript variants encoding different isoforms have been described for this gene.
Applications:	ELISA, IHC
Name of antibody:	CASP10
Immunogen:	Synthetic peptide of human CASP10
Full name:	caspase 10
Synonyms:	MCH4; ALPS2; FLICE2
SwissProt:	Q92851
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer; Human ovarian cancer
IHC Recommend dilution:	20-100



