

## 兔抗 DIRAS1 多克隆抗体

中文名称：兔抗 DIRAS1 多克隆抗体

英文名称：Anti-DIRAS1 rabbit polyclonal antibody

别名：RIG; GBTS1; Di-Ras1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DIRAS1

反应种属：Human


标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	DIRAS1 (DIRAS family, GTP-binding Ras-like 1), also known as distinct subgroup of the Ras family member 1, RIG (Ras-related inhibitor of cell growth), small GTP-binding tumor suppressor 1 or GBTS1, is a 198 amino acid cell membrane protein expressed at high levels in brain and heart. DIRAS1 displays low GTPase activity and is encoded by a gene that maps to human chromosome 19p13.3.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DIRAS1
<b>Immunogen:</b>	Synthetic peptide of human DIRAS1

<b>Full name:</b>	DIRAS family, GTP-binding RAS-like 1
<b>Synonyms :</b>	RIG; GBTS1; Di-Ras1
<b>SwissProt:</b>	O95057
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human cervical cancer and Human esophagus cancer
<b>IHC Recommend dilution:</b>	50-200



**mlbio** 梅联生物  
Good elisakit producers