

兔抗 DIRAS3 多克隆抗体

中文名称：兔抗 DIRAS3 多克隆抗体

英文名称： Anti-DIRAS3 rabbit polyclonal antibody

别名： DIRAS family GTPase 3; ARHI; NOEY2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： DIRAS3

反应种属： Human, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the ras superfamily. This gene is imprinted gene with monoallelic expression of the paternal allele which is associated with growth suppression. The encoded protein acts as a tumor suppressor whose function is abrogated in many ovarian and breast cancers. This protein may also play a role autophagy in certain cancer cells by regulating the autophagosome initiation complex.
Applications:	ELISA, WB
Name of antibody:	DIRAS3
Immunogen:	Fusion protein of human DIRAS3

Full name:	DIRAS family GTPase 3
Synonyms:	ARHI; NOEY2
SwissProt:	O95661
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	26 kDa
WB Positive control:	Rat brain tissue and Rat liver tissue lysates
WB Recommended dilution:	200-1000

