

兔抗 CASQ1 多克隆抗体

中文名称：兔抗 CASQ1 多克隆抗体

英文名称：Anti-CASQ1 rabbit polyclonal antibody

别名：calsequestrin 1; CASQ; PDIB1; VMCQA

相关类别：一抗

抗原：CASQ1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes the skeletal muscle specific member of the calsequestrin protein family. Calsequestrin functions as a luminal sarcoplasmic reticulum calcium sensor in both cardiac and skeletal muscle cells. This protein, also known as calmitine, functions as a calcium regulator in the mitochondria of skeletal muscle. This protein is absent in patients with Duchenne and Becker types of muscular dystrophy.
Applications:	ELISA, IHC
Name of antibody:	CASQ1
Immunogen:	Fusion protein of human CASQ1

Full name:	calsequestrin 1
Synonyms:	CASQ; PDIB1; VMCQA
SwissProt:	P31415
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-300

