

兔抗 DIXDC1 多克隆抗体

中文名称：兔抗 DIXDC1 多克隆抗体

英文名称： Anti-DIXDC1 rabbit polyclonal antibody

别名： C3orf41

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： DIXDC1

宿主： Rabbit

反应种属： Human, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a positive regulator of the Wnt signaling pathway. The encoded protein is found associated with gamma tubulin at the centrosome. Alternative splicing results in multiple transcript variants encoding different isoforms.
Applications:	ELISA, IHC
Name of antibody:	DIXDC1
Immunogen:	Synthetic peptide corresponding to a region derived from 253-268 amino acids of human DIXDC1
Full name:	DIX domain containing 1
Synonyms:	CCD1

SwissProt:	Q155Q3
ELISA Recommended dilution:	5000-10000
IHC Positive control:	Human brain;
IHC Recommended dilution:	10-50

