

兔抗 DKK1 多克隆抗体

中文名称：兔抗 DKK1 多克隆抗体

英文名称：Anti-DKK1 rabbit polyclonal antibody

别名：SK; DKK-1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DKK1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein that is a member of the dickkopf family. It is a secreted protein with two cysteine rich regions and is involved in embryonic development through its inhibition of the WNT signaling pathway. Elevated levels of DKK1 in bone marrow plasma and peripheral blood is associated with the presence of osteolytic bone lesions in patients with multiple myeloma.
Applications:	ELISA, IHC
Name of antibody:	DKK1
Immunogen:	Synthetic peptide of human DKK1
Full name:	dickkopf WNT signaling pathway inhibitor 1

Synonyms :	SK; DKK-1
SwissProt:	O94907
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human brain
IHC Recommend dilution:	25-100

