

## 兔抗 DKK3 多克隆抗体

中文名称：兔抗 DKK3 多克隆抗体

英文名称： Anti-DKK3 rabbit polyclonal antibody

别名： RIG; REIC

抗原： DKK3

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein that is a member of the dickkopf family. The secreted protein contains two cy steine rich regions and is involved in embryonic deve lopment through its interactions with the Wnt signali ng pathway. The expression of this gene is decreased in a variety of cancer cell lines and it may function a s a tumor suppressor gene. Alternative splicing result s in multiple transcript variants encoding the same pr otein.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DKK3

<b>Immunogen:</b>	Fusion protein of human DKK3
<b>Full name:</b>	dickkopf WNT signaling pathway inhibitor 3
<b>Synonyms :</b>	RIG; REIC
<b>SwissProt:</b>	Q9UBP4
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human thyroid cancer and human gastric cancer
<b>IHC Recommend dilution:</b>	25-100

