

兔抗 CASZ1 多克隆抗体

中文名称: 兔抗 CASZ1 多克隆抗体

英文名称: Anti-CASZ1 rabbit polyclonal antibody

别名: CST; SRG; CAS11; ZNF693; dJ734G22.1

储存: 冷冻 (-20℃)

抗原: CASZ1

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a zinc finger transcription factor. The encoded protein may function as a tumor suppressor, and single nucleotide polymorphisms in this gene are associated with blood pressure variation. Alternative splicing results in multiple transcript variants that encode different protein isoforms.
Applications:	ELISA, IHC
Name of antibody:	CASZ1
Immunogen:	Synthetic peptide of human CASZ1
Full name:	castor zinc finger 1

Synonyms :	CST; SRG; CAS11; ZNF693; dJ734G22.1
SwissProt:	Q86V15
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	25-100

