

兔抗 DKKL1 多克隆抗体

中文名称： 兔抗 DKKL1 多克隆抗体

英文名称： Anti-DKKL1 rabbit polyclonal antibody

别名： SGY; CT34; SGY1; SGY-1

相关类别： 一抗

抗原： DKKL1

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The dickkopf protein family interacts with the Wnt signaling pathway and its members are characterized by two conserved cysteine-rich domains. This gene encodes a secreted protein that has low sequence similarity to the dickkopf-3 protein. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	DKKL1
Immunogen:	Synthetic peptide corresponding to a region derived from 202-214 amino acids of human DKKL1

Full name:	dickkopf like acrosomal protein 1
Synonyms:	SGY; CT34; SGY1; SGY-1
SwissProt:	Q9UK85
ELISA Recommended dilution:	5000-10000
IHC Positive control:	Human esophagus cancer;Human liver cancer
IHC Recommended dilution:	25-100



