

## 兔抗 CATSPER1 多克隆抗体

中文名称：兔抗 CATSPER1 多克隆抗体

英文名称： Anti-CATSPER1 rabbit polyclonal antibody

别名： SPGF7; CATSPER

相关类别： 一抗

抗原： CATSPER1

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Calcium ions play a primary role in the regulation of sperm motility. This gene belongs to a family of putative cation channels that are specific to spermatozoa and localize to the flagellum. The protein family features a single repeat with six membrane-spanning segments and a predicted calcium-selective pore region.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CATSPER1
<b>Immunogen:</b>	Fusion protein of human CATSPER1
<b>Full name:</b>	cation channel, sperm associated 1
<b>Synonyms :</b>	SPGF7; CATSPER

<b>SwissProt:</b>	Q8NEC5
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human thyroid cancer and Human prostate cancer
<b>IHC Recommend dilution:</b>	25-100

