

## 兔抗 CBX6 多克隆抗体

中文名称：兔抗 CBX6 多克隆抗体

英文名称： Anti-CBX6 rabbit polyclonal antibody

储 存： 冷冻（-20℃）

抗 原： CBX6

宿 主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	CBX6 (chromobox protein homolog 6) is a 412 amino acid protein that is a component of the chromatin-associated PcG complex. Through interaction with PcG, CBX6 modifies chromatin and maintains the transcriptionally repressive state of genes. Members of the chromobox domain protein family are characterized by the presence of a 37 amino acid chromobox (CHRromatin Organization MODifier) domain. This domain plays a mechanistic role in targeting chromodomain proteins to specific regions of the nucleus.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	CBX6
<b>Immunogen:</b>	Fusion protein of human CBX6
<b>Full name:</b>	chromobox homolog 6

<b>SwissProt:</b>	O95503
<b>ELISA Recommended dilution:</b>	2000-5000
<b>WB Predicted band size:</b>	44 kDa
<b>WB Positive control:</b>	Hela and 293T cells
<b>WB Recommended dilution:</b>	500-2000

