

兔抗 CCAR1 多克隆抗体

中文名称：兔抗 CCAR1 多克隆抗体

英文名称： Anti-CCAR1 rabbit polyclonal antibody

别名： RP11-437A18.1

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： CCAR1

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Cell division cycle and apoptosis regulator protein 1 is a protein that in humans is encoded by the CCAR1 gene. The deduced 1,150-amino acid CARP1 protein contains a putative cold-shock protein domain expected to bind RNA and a putative DNA-binding motif predicted to be involved in chromosomal organization. Immunohistochemical analysis detected CARP1 in the perinuclear compartment of human breast cancer cells, and the protein had an apparent molecular mass of 130 kD by SDS-PAGE.
Applications:	ELISA, IHC

Name of antibody:	CCAR1
Immunogen:	Fusion protein of human CCAR1
Full name:	cell division cycle and apoptosis regulator 1
Synonyms:	RP11-437A18.1
SwissProt:	Q8IX12
IHC positive control:	Human gastric cancer and human liver cancer
IHC Recommend dilution:	50-200

