

## 兔抗 DLL3 多克隆抗体

中文名称：兔抗 DLL3 多克隆抗体

英文名称：Anti-DLL3 rabbit polyclonal antibody

别名：SCDO1

相关类别：一抗

储存：冷冻（-20℃）

抗原：DLL3

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the delta protein ligand family. This family functions as Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain. Mutations in this gene cause autosomal recessive spondylocostal dysostosis 1. Two transcript variants encoding distinct isoforms have been identified for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DLL3
<b>Immunogen:</b>	Synthetic peptide of human DLL3
<b>Full name:</b>	delta like canonical Notch ligand 3

<b>Synonyms:</b>	SCDO1
<b>SwissProt:</b>	Q9NYJ7
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer;Human prostate cancer
<b>IHC Recommend dilution:</b>	20-100



