

兔抗 DLL4 多克隆抗体

中文名称：兔抗 DLL4 多克隆抗体

英文名称： Anti-DLL4 rabbit polyclonal antibody

别名： hdelta2

相关类别： 一抗

抗原： DLL4

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is a homolog of the Drosophila delta gene. The delta gene family encodes Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain.
Applications:	ELISA, IHC
Name of antibody:	DLL4
Immunogen:	Fusion protein of human DLL4
Full name:	delta like canonical Notch ligand 4
Synonyms:	AOS6; delta4; hdelta2
SwissProt:	Q9NR61
ELISA Recommended dilution:	2000-5000

IHC positive control:	Human breast cancer and Human thyroid cancer
IHC Recommend dilution:	50-200

