

## 兔抗 CCDC12 多克隆抗体

中文名称：兔抗 CCDC12 多克隆抗体

英文名称：Anti-CCDC12 rabbit polyclonal antibody

别名：coiled-coil domain containing 12;

相关类别：一抗

抗原：CCDC12

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

**Background:**

CCDC12, also known as FLJ39430, FLJ40801 or MGC23918, is a 166 amino acid protein encoded by a gene mapping to human chromosome 3. Chromosome 3 is made up of about 214 million bases encoding over 1,100 genes. Notably, there is a chemokine receptor gene cluster and a variety of human cancer related loci on chromosome 3. Particular regions of the chromosome 3 short arm are deleted in many types of cancer cells. Key tumor suppressing genes on chromosome 3 encode apoptosis mediator RASSF1, cell migration regulator HYAL1 and angiogenesis suppressor S

	EMA3B. Marfan Syndrome, porphyria, von Hippel-Lindau syndrome, osteogenesis imperfecta and Charcot-Marie-Tooth Disease are a few of the numerous genetic diseases associated with chromosome 3.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CCDC12
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	coiled-coil domain containing 12
<b>SwissProt:</b>	Q8WUD4
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer and human ovarian cancer
<b>IHC Recommend dilution:</b>	50-300



