

兔抗 CCDC47 多克隆抗体

中文名称：兔抗 CCDC47 多克隆抗体

英文名称： Anti-CCDC47 rabbit polyclonal antibody

别名： coiled-coil domain containing 47; GK001; MSTP041

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CCDC47

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The coiled-coil domain is a common protein motif that is often involved in protein oligomerization and is found in proteins such as transcription factors and intermediate filaments. The CCDC47 gene maps to chromosome 17 at 17q23.3. Little is known about this single-pass membrane protein except that the coiled-coil domain is within the cytosolic domain near the carboxy terminus.
Applications:	ELISA, WB, IHC
Name of antibody:	CCDC47
Immunogen:	Fusion protein of human CCDC47

Full name:	coiled-coil domain containing 47
Synonyms:	GK001; MSTP041
SwissProt:	Q96A33
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer and Human tonsil
IHC Recommend dilution:	30-150
WB Predicted band size:	56 kDa
WB Positive control:	Mouse heart tissue, Mouse liver tissue, NIH/3T3 cell, Human bladder transitional cell carcinoma grade 2-3 tissue, A172 cell, Human placenta tissue, K562 cell lysates
WB Recommended dilution:	500-2000



