

## 兔抗 DLX5 多克隆抗体

中文名称：兔抗 DLX5 多克隆抗体

英文名称： Anti-DLX5 rabbit polyclonal antibody

别名： SHFM1D

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： DLX5

宿主： Rabbit

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes a member of a homeobox transcription factor gene family similar to the Drosophila distal-less gene. The encoded protein may play a role in bone development and fracture healing. Mutation in this gene, which is located in a tail-to-tail configuration with another member of the family on the long arm of chromosome 7, may be associated with split-hand/split-foot malformation.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	DLX5
<b>Immunogen:</b>	Synthesized peptide derived from internal of human

	DLX5.
<b>Full name:</b>	distal-less homeobox 5
<b>Synonyms :</b>	SHFM1D
<b>SwissProt:</b>	P56178
<b>IHC positive control:</b>	Human lung carcinoma tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	32 kDa
<b>WB Positive control:</b>	Hela cells, 293 cells and COLO205 cells lysates
<b>WB Recommended dilution:</b>	500-3000



