

兔抗 CCDC102A 多克隆抗体

中文名称：兔抗 CCDC102A 多克隆抗体

英文名称：Anti-CCDC102A rabbit polyclonal antibody

别名：coiled-coil domain containing 102A;

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CCDC102A

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

| | |
|------------------------------------|---|
| Background: | The function of this protein remains unknown. |
| Applications: | ELISA, IHC |
| Name of antibody: | CCDC102A |
| Immunogen: | Synthetic peptide of human CCDC102A |
| Full name: | coiled-coil domain containing 102A |
| SwissProt: | Q96A19 |
| ELISA Recommended dilution: | 5000-10000 |
| IHC positive control: | Human liver cancer |
| IHC Recommend dilution: | 30-150 |

