

兔抗 DMC1 多克隆抗体

中文名称：兔抗 DMC1 多克隆抗体

英文名称： Anti-DMC1 rabbit polyclonal antibody

别名： DMC1H; LIM15; dJ199H16.1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： DMC1

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the superfamily of recombinases (also called DNA strand-exchange proteins). Recombinases are important for repairing double-strand DNA breaks during mitosis and meiosis. This protein, which is evolutionarily conserved, is reported to be essential for meiotic homologous recombination and may thus play an important role in generating diversity of genetic information. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	DMC1

Immunogen:	Full length fusion protein(165FULL)
Full name:	DNA meiotic recombinase 1
Synonyms:	DMC1H; LIM15; dJ199H16.1
SwissProt:	Q14565
ELISA Recommended dilution:	5000-10000
IHC Positive control:	Human gastric cancer;Human lung cancer
IHC Recommended dilution:	30-150



