

## 兔抗 CCDC113 多克隆抗体

中文名称：兔抗 CCDC113 多克隆抗体

英文名称：Anti-CCDC113 rabbit polyclonal antibody

别名：coiled-coil domain containing 113; HSPC065

相关类别：一抗

抗原：CCDC113

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Component of centriolar satellites contributing to primary cilium formation.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CCDC113
<b>Immunogen:</b>	Fusion protein of human CCDC113
<b>Full name:</b>	coiled-coil domain containing 113
<b>Synonyms:</b>	HSPC065
<b>SwissProt:</b>	Q9H0I3
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human ovarian cancer and human esophagus cancer

IHC Recommend dilution:

30-150



