

兔抗 DMRTA2 多克隆抗体

中文名称：兔抗 DMRTA2 多克隆抗体

英文名称：Anti-DMRTA2 rabbit polyclonal antibody

别名：DMRT like family A2

抗原：DMRTA2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May be involved in sexual development.
Applications:	ELISA, IHC
Name of antibody:	DMRTA2
Immunogen:	Synthetic peptide of human DMRTA2
Full name:	DMRT like family A2
SwissProt:	Q96SC8
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human esophagus cancer
IHC Recommend dilution:	50-200

