

## 兔抗 CCIN 多克隆抗体

中文名称：兔抗 CCIN 多克隆抗体

英文名称：Anti-CCIN rabbit polyclonal antibody

别名：calicin; BTBD20; KBTBD14

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CCIN

反应种属：Human, Mouse, Rat

标记物：Unconjugate 克

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a basic protein of the sperm head cytoskeleton. This protein contains kelch repeats and a BTB/POZ domain and is necessary for normal morphology during sperm differentiation. This gene is intronless.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CCIN
<b>Immunogen:</b>	Synthetic peptide of human CCIN
<b>Full name:</b>	calicin
<b>Synonyms:</b>	BTBD20; KBTBD14
<b>SwissProt:</b>	Q13939

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	10-50

