

## 兔抗 DNA2 多克隆抗体

中文名称：兔抗 DNA2 多克隆抗体

英文名称： Anti-DNA2 rabbit polyclonal antibody

别名： DNA replication helicase/nuclease 2; DNA2L; hDNA2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： DNA2

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the DNA2/NAM7 helicase family. The encoded protein is a conserved helicase/nuclease involved in the maintenance of mitochondrial and nuclear DNA stability. Mutations in this gene are associated with autosomal dominant progressive external ophthalmoplegia-6 (PEOA6) and Seckel syndrome 8. Alternatively spliced transcript variants have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DNA2
<b>Immunogen:</b>	Fusion protein of human DNA2
<b>Full name:</b>	DNA replication helicase/nuclease 2

<b>Synonyms:</b>	DNA2L; hDNA2
<b>SwissProt:</b>	P51530
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human breast cancer
<b>IHC Recommend dilution:</b>	50-100

