

兔抗 DNaAF5 多克隆抗体

中文名称：兔抗 DNaAF5 多克隆抗体

英文名称：Anti-DNaAF5 rabbit polyclonal antibody

别名：CILD18; HEATR2

相关类别：一抗

抗原：DNaAF5

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is essential for the preassembly or stability of axonemal dynein arms, and is found only in organisms with motile cilia and flagella. Mutations in this gene are associated with primary ciliary dyskinesia-18, a disorder characterized by abnormalities of motile cilia. Alternatively spliced transcript variants have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	DNaAF5
Immunogen:	Fusion protein of human DNaAF5
Full name:	dynein, axonemal, assembly factor 5

Synonyms :	CILD18; HEATR2
SwissProt:	Q86Y56
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	25-100

