

## 兔抗 CCL2 多克隆抗体

中文名称：兔抗 CCL2 多克隆抗体

英文名称： Anti-CCL2 rabbit polyclonal antibody

别名： HC11, MCAF, MCP1, MCP-1, SCYA2, GDCF-2, SMC-CF, HSMCR30

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： CCL2

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

This gene is one of several cytokine genes clustered on the q-arm of chromosome 17. Cytokines are a family of secreted proteins involved in immunoregulatory and inflammatory processes. The protein encoded by this gene is structurally related to the CXC subfamily of cytokines. Members of this subfamily are characterized by two cysteines separated by a single amino acid. This cytokine displays chemotactic activity for monocytes and basophils but not for neutrophils or eosinophils. It has been implicated in the pathogenesis of diseases characterized by monocytic infiltrates, like psoriasis, rhe

	umatoid arthritis and atherosclerosis. It binds to chemokine receptors CCR2 and CCR4.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CCL2
<b>Immunogen:</b>	Synthetic peptide of human CCL2
<b>Full name:</b>	C-C motif chemokine 2
<b>Synonyms :</b>	HC11, MCAF, MCP1, MCP-1, SCYA2, GDCF-2, SMC-CF, HSMCR30
<b>SwissProt:</b>	P13500
<b>ELISA Recommended dilution:</b>	500-5000
<b>IHC positive control:</b>	Human breast cancer
<b>IHC Recommend dilution:</b>	25-50

